# UMA PHARMATECH MACHINERY THE FUTURE OF BIOTECH



"Designed To Work Efficiently"

"Minimum power consumption"

"Energy Efficient Impellers"

"Customized Mixing Equipment" Address : 61/62, Yamuna Industrial Park, Vanch Dhamatvan Road, Dhamatvan, Pin:382435, Ahmedabad | Gujarat | India Contact: 9033404554/9726923885/8866137364 Email:umapharmatechmachinery@gmail.com

# UMA PHARMATECH MACHINERY THE FUTURE OF BIOTECH



## Who We Are?

We Are On Of The Leading Manufacturer for Biotech and Life Science Equipment's in Gujarat. We Have More Then 10+ Years Of Experience in Design Of Machinery and Plant Setup. We Have Well Experienced Design Team For Manufacturing Of Various Lab Equipment's. We Provide Turnkey Solution For Bio fertilizer, Bio pesticides, Bio plastics, Probiotics, API With Utilities and Automation.

# Why choose us?

Designated to Work Efficiently without Compromising in Aseptic Fermentation Technology under one roof.

In-house Manufacturing Facility Like EOT Crane(10 TON), Plasma Cutting Machine, Argon Welding Machinery, X-ray Welding Testing, Hydro Testing, Accurate Machining, Drilling any Many More.

#### UMA PHARMATECH MACHINERY

The Feature Of Biotech

M: +91-9726923885 Email: info@ssfermenterbioreactor.com Web: www.umapharmatechmachinery.com Address: Plot No. 61, 62, Yamuna Industrial Park, Dhamatvan, Ahmedabad-382435, Gujarat, India

## **Our Product**

- ➡ Lab Fermenter(5 Litre To 20 Litre)
- Pilot scale fermenter((50 Litre To 1000)
- production Scale Fermenter(1.5kl to 60kl)
- ➡ Vertical Double Wall Autoclave Sterilizer
- Laboratory and Heavy Duty Rotary Shaker
- BOD Incubator
- SS Reactor and SS Tanks
- ➡ Powder Mixer Machine
- Turnkey Plant For API, Bio fertilizer, Probiotics

## Utilities

- Fermenter(Insitu/ Autoclavable)
- Oil Free Air Compressor
- Electric Steam Generator/ Vertical Autoclave
- 🔷 Chillier Unit
- Electric Control Panel
  - SS Skid For Line Distribution With Utility For Insitu Fermenter

### **Features**

- Lab Fermenter(5 Litre To 20 Litre)
  - Autoclavable and Insitu Model Available
  - Autoclavable and Insitu Model Available
- Automatic Temperature Controlling System
- Automatic Sterilization System
- Automatic Ph Controlling System
  - PLC Based Or Microprocessor Based Control
- Agitation Speed Controlling Knob

# UMA PHARMATECH MACHINERY

The Feature Of Biotech



### Standard Model of rotary shaker/ Orbital Shaker

#### Heavy Duty Industrial Model

Platform Size(mm) 1250 x 1250							
Flask Capacity(ml) 500 1000 2000 5000							
Number Of Flask 81 64 36 25							

Platform Size(mm) 1000 x 1000								
Flask Capacity(ml) 500 1000 2000 5000								
Number Of Flask4936259								

#### Laboratory scale Model

Platform Size(mm) 850 x 850						
Flask Capacity(ml) 250 500 1000 2000 5000						
Number Of Flask 49 36 25 16 9						

Platform Size(mm) 600 x 600						
Flask Capacity(ml) 250 500 1000 2000 5000						
Number Of Flask	25	16	9	9	4	

### UMA PHARMATECH MACHINERY

The Feature Of Biotech

### Standard Model of rotary shaker/ Orbital Shaker

Total	Inner Chamber Size(cm)	Heater				
Volume(litre)		Capacity(Watt)				
Single Wa	II Vertical Autoclave Without Foot L	ifting Facility				
21	30 dia x 30 height	2000				
35	30 dia x 50 height	2000				
53	37 dia x 50 height	3000				
75	37 dia x 60 height	3000				
Double \	Double Wall Vertical Autoclave With Foot Lifting Facility					
100	47 dia x 58 height	6000				
150	47 dia x 87 height	9000				
210	55 dia x 90 height	12000				

## UMA PHARMATECH MACHINERY

The Feature Of Biotech



#### **BOD INCUBATOR WITH HEATING AND COOLING**

#### LAB BOD INCUBATOR

Working Size( H x W x D)(cm)	Liters	Trays
45 x 45 x 45	92	2
60 x 45 x 45	121	2
70 x 50 x 50	175	2
60 x 60 x 60	216	3
90 x 60 x 60	324	3

UMA PHARMATECH MACHINERY

**The Feature Of Biotech** 



#### **Electric Seam Boiler/ Steam Generator**

#### Industrial Model

Model No	Steam Production (kg/hr)	Max. Design pressure	Max. Operating pressure	Connected Load(Kw)
		(Bar)	(Bar)	
UPMB 12	15	6	3	12
UPMB 15	19	6	3	15
UPMB 18	23	6	3	18
UPMB 24	31	6	3	24
UPMB 27	35	6	3	27
UPMB 30	39	6	3	30
UPMB 36	47	6	3	36
UPMB 48	62	6	3	48
UPMB 60	78	6	3	60

UMA PHARMATECH MACHINERY

**The Feature Of Biotech** 

M: +91-9726923885 Email: info@ssfermenterbioreactor.com Web: www.umapharmatechmachinery.com Address: Plot No. 61, 62, Yamuna Industrial Park, Dhamatvan, Ahmedabad-382435, Gujarat, India

#### Features

- Automatic Temperature Controlling and Pressure Controlling
- Industrial Heavy Duty Immersion Heating Element
- Automatic Feed Water Available On Demand
- Glass Level Indicator For Water Level
- Inner Shell Made Of ss304/Ms On Demand
- Outer Body and Cover Made Of Mild Steel Painted





#### **Ribbon Blender**

Model	Total	Working	Working	Output	Motor(hp
	Volume	Volume(LTR)	Volume(	RPM	) Approx.
	(LTR)		kg)		
UPMR50	55	40	50	50-60	1/2
UPMR100	105	85	100	40-50	2/3
UPMR200	215	170	200	35-45	2/3
UPMR300	308	245	300	30-45	3/5
UPMR500	517	415	500	30-45	5/7.5
UPMR1050	1060	845	1050	25-40	5/7.5/10

Working Volume (Kg) calculated based on Density Og Material 1.25 g/cm<sup>3</sup>

#### UMA PHARMATECH MACHINERY

**The Feature Of Biotech** 

QUALITY IS Our commitment

# **Our Customer**



















# GreenMax AgroTech

Coimbatore, Tamil Nadu









